


GUIDEZILLA™ II

Guide Extension Catheter




General Specifications

Indication:	The GUIDEZILLA II Guide Extension Catheter is intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.
Available Sizes:	6F, 7F, 8F, and 6F Long*
Coating:	Z-Glide™ hydrophilic coating on distal guide segment
Distal Guide Length:	6F*, 7F*, 8F* = 25 cm 6F Long* = 40 cm
Proximal Shaft:	Stainless steel hypotube
Collar Type:	Platinum iridium collar embedded in polymer

*6F = 2 mm, 7F = 2.33 mm, 8F = 2.67 mm

GUIDEZILLA™ II

Guide Extension Catheter



GUIDEZILLA II Guide Extension Catheter Ordering Information

Size	GTIN	Ref/Catalog Number	Compatible Guide Catheter	Inner Diameter	Outer Diameter
6F*	08714729939450	H749 393351506 0	6F I.D. \geq 0.070" (1.78 mm)	0.057" (1.45 mm)	0.067" (1.71 mm)
7F*	08714729939474	H749 393351507 0	7F I.D. \geq 0.078" (1.98 mm)	0.063" (1.60 mm)	0.073" (1.86 mm)
8F*	08714729939481	H749 393351508 0	8F I.D. \geq 0.088" (2.24 mm)	0.072" (1.83 mm)	0.083" (2.11 mm)
6F Long*	08714729939467	H749 3933515061 0	6F I.D. \geq 0.070" (1.78 mm)	0.057" (1.45 mm)	0.067" (1.71 mm)

*6F = 2 mm, 7F = 2.33 mm, 8F = 2.67 mm

Boston Scientific
Advancing science for life™

Caution: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for use or distribution outside the US, France, and Japan only. All trademarks are the property of their respective owners.

IC-455611-AA-EU MAR2017

www.bostonscientific.eu

© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
DINCAR2714EA