


Difficult Terrain. TERUMO Territory.™



A broad range of
guiding sheath shapes
and sizes to meet your
interventional and
diagnostic challenges

Pinnacle® Destination®
Guiding Sheaths

Terumo Pinnacle® Destination® Guiding Sheaths safe

- **Enhances visualization for precise positioning**

- Triple radiopacity (sheath, dilator, and gold coil)

- **Maintains shape and positioning**

- Outstanding durability, even in the longest procedures

- **Provides lesion access**

- Exceptional trackability, even in the most challenging anatomy


- **Accommodates anatomical variations in vessel access**

- Assorted tip shapes

ely and efficiently get you there and keep you there.

- **Allows smooth transition and easy penetration**

- Unique atraumatic tip minimizes potential for vessel damage
- Smooth guidewire-to-dilator and dilator-to-sheath transitions



- **Resists kinking**

- Stainless steel coil for consistent reliability

- **Provides smooth movement**
















































- Coil reinforced tubing and a PTFE inner layer that minimizes friction

From the creators of the Glidewire®

**Because there's no such thing as a "standard" procedure,
you're always in TERUMO Territory.**

Pinnacle Destination Guiding Sheaths

Packed 1 per box and includes dilator. All dilators are 0.038" wire compatible.

French Size	Product Code	Length (cm)	Hydrophilic Coating Length (cm)	Valve Type	Curve Style	Curve Shape
5 Fr	54-54502	45	5	Cross-Cut	Hockey Stick	
	54-54505	45	5	Cross-Cut	LIMA	
	54-54503	45	5	Cross-Cut	Multi-Purpose	
	54-54504	45	5	Cross-Cut	RDC	
	54-54501	45	35	Cross-Cut	Straight	
6 Fr	RSR02	45	5	Cross-Cut	Hockey Stick	
	RSR14	45	5	Cross-Cut	LIMA	
	RSR03	45	5	Cross-Cut	Multi-Purpose	
	RSR13	45	5	Cross-Cut	RDC	
	RSR01	45	5	Cross-Cut	Straight	
	RSR08	45	5	Tuohy-Borst	Hockey Stick	
	RSR18	45	5	Tuohy-Borst	LIMA	
	RSR09	45	5	Tuohy-Borst	Multi-Purpose	
	RSR17	45	5	Tuohy-Borst	RDC	
	RSR07	45	5	Tuohy-Borst	Straight	
	54-64501	45	35	Cross-Cut	Straight	
	54-64506	45	35	Tuohy-Borst	Straight	
	RSP01	65	35	Cross-Cut	Straight	
	RSP03	65	35	Tuohy-Borst	Straight	
	RSC07	90	15	Cross-Cut	Multi-Purpose	
	RSC05	90	15	Cross-Cut	Straight	
	RSC03	90	15	Tuohy-Borst	Multi-Purpose	
	RSC01	90	15	Tuohy-Borst	Straight	
7 Fr	RSR05	45	5	Cross-Cut	Hockey Stick	
	RSR16	45	5	Cross-Cut	LIMA	
	RSR06	45	5	Cross-Cut	Multi-Purpose	
	RSR15	45	5	Cross-Cut	RDC	
	RSR04	45	5	Cross-Cut	Straight	
	RSR11	45	5	Tuohy-Borst	Hockey Stick	
	RSR20	45	5	Tuohy-Borst	LIMA	
	RSR12	45	5	Tuohy-Borst	Multi-Purpose	
	RSR19	45	5	Tuohy-Borst	RDC	
	RSR10	45	5	Tuohy-Borst	Straight	
	54-74501	45	35	Cross-Cut	Straight	
	54-74506	45	35	Tuohy-Borst	Straight	
	RSP02	65	35	Cross-Cut	Straight	
	RSP04	65	35	Tuohy-Borst	Straight	
	RSC08	90	15	Cross-Cut	Multi-Purpose	
	RSC06	90	15	Cross-Cut	Straight	
	RSC04	90	15	Tuohy-Borst	Multi-Purpose	
	RSC02	90	15	Tuohy-Borst	Straight	
8 Fr	54-84501	45	35	Cross-Cut	Straight	
	54-84506	45	35	Tuohy-Borst	Straight	
	54-86501	65	35	Cross-Cut	Straight	
	54-86506	65	35	Tuohy-Borst	Straight	
	54-89001	90	60	Cross-Cut	Straight	
	54-89006	90	60	Tuohy-Borst	Straight	

For more information, call 800.862.4143 to speak to an Inside Sales Customer Care Representative, or visit www.terumo.com
Terumo Interventional Systems • 2101 Cottontail Lane • Somerset, NJ 08873 • Fax: 800.411.5870

©2009 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit #TIS122r1 07/09