


AVIATOR[™] Plus

PTA Dilatation Catheter
– for speed and
efficiency in carotid
and renal peripheral
endovascular
interventions



AVIATOR™ Plus

The fully dedicated 0.014" rapid-exchange



Outstanding deliverability and crossability

- Ultra low profile balloon with new optimal pleating and folding

Excellent pushability


- Shaft designed to overcome the most challenging anatomies

Superb trackability

- Smooth tip-to-wire transition

Remarkable pressure resistance


- Cordis DURALYN™ balloon technology providing
 - Excellent puncture resistance
 - Short shoulders
 - Controlled compliance



balloon for carotid and renal arteries

Designed for ease of use and speed combined in a low profile, efficient system


Balloon Deflation Times*



* Cordis data on file.
** Viatrac is a registered trademark of Abbott
*** Sterling OTW is a registered trademark of Boston Scientific
****Submarine is a registered trademark of Invatec

Product Specifications

- Balloon diameter range: 4 - 7 mm
- Rated Burst Pressure: up to 14 atm
- DURALYN™ balloon material
- Shaft diameter: 3.3 F
- Usable catheter length: 142 cm
- Compatible with 0.014" interventional guidewire
- Ideal for carotid post-dilatation and renal procedures



AVIATOR™ Plus

PTA Dilatation Catheter – for speed and efficiency in carotid and renal peripheral endovascular interventions

Ordering Information

Balloon OD x Length (mm x cm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Guiding Catheter Fit (F)	Catalogue Number
4 x 1.5	10	14	4	6	424-4015 W
4 x 2	10	14	4	6	424-4020 W
4 x 3	10	14	4	6	424-4030 W
4 x 4	10	14	4	6	424-4040 W
4.5 x 1.5	10	14	4	6	424-4515 W
4.5 x 2	10	14	4	6	424-4520 W
4.5 x 3	10	14	4	6	424-4530 W
4.5 x 4	10	14	4	6	424-4540 W
5 x 1.5	10	14	4	6	424-5015 W
5 x 2	10	14	4	6	424-5020 W
5 x 3	10	14	4	6	424-5030 W
5 x 4	10	14	4	6	424-5040 W
5.5 x 1.5	10	14	4	6	424-5515 W
5.5 x 2	10	14	4	6	424-5520 W
5.5 x 3	10	14	4	6	424-5530 W
5.5 x 4	10	14	4	6	424-5540 W
6 x 1.5	10	14	4	6	424-6015 W
6 x 2	10	14	4	6	424-6020 W
6 x 3	10	14	4	6	424-6030 W
6 x 4	10	14	4	6	424-6040 W
7 x 1.5	10	12	4	6	424-7015 W
7 x 2	10	12	4	6	424-7020 W
7 x 3	10	12	4	6	424-7030 W
7 x 4	10	12	4	7	424-7040 W

Important Information: Prior to use, refer to the "instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

© Cordis Europe, a division of Johnson & Johnson Medical N.V. • May 2012 • 2E-800-1655-3 • Tel: +32 2 746 30 00 • Fax: +32 2 746 34 05