


We Are Helping You Save Patients' Limbs by Expanding the Portfolio


INTRODUCING

Cordis POWERFLEX[®] Pro

0.035" PTA Dilatation Catheter



**Advanced Crossing Ability.
Remarkable Versatility.**




Cordis **POWERFLEX[®] Pro**

0.035" PTA Dilatation Catheter



Profile, Power, and Performance for Treating Diffuse and Challenging Lower Extremity Lesions.

Sizes up to 220mm allow for the treatment of long and diffuse lesions with single dilatation.



A New Addition to the Lower Extremity Solutions Portfolio

Ordering Information

80cm Shaft Length

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	CSI Fit F
4400302S	3	20	5
4400304S	3	40	5
4400306S	3	60	5
4400308S	3	80	5
4400310S	3	100	5
4400312S	3	120	5
4400315S	3	150	5
4400322S	3	220	5
4400402S	4	20	5
4400404S	4	40	5
4400406S	4	60	5
4400408S	4	80	5
4400410S	4	100	5
4400412S	4	120	5
4400415S	4	150	5
4400422S	4	220	5
4400502S	5	20	5
4400503S	5	30	5
4400504S	5	40	5
4400506S	5	60	5
4400508S	5	80	5
4400510S	5	100	5
4400512S	5	120	5
4400515S	5	150	5
4400522S	5	220	5
4400602S	6	20	5
4400603S	6	30	5
4400604S	6	40	5
4400606S	6	60	5
4400608S	6	80	5
4400610S	6	100	5

80cm Shaft Length Continued

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	CSI Fit F
4400612S	6	120	5
4400615S	6	150	5
4400622S	6	220	5
4400702S	7	20	5
4400703S	7	30	5
4400704S	7	40	5
4400706S	7	60	5
4400708S	7	80	5
4400710S	7	100	6
4400802S	8	20	5
4400803S	8	30	5
4400804S	8	40	5
4400806S	8	60	5
4400808S	8	80	5
4400810S	8	100	6
4400902S	9	20	6
4400903S	9	30	6
4400904S	9	40	6
4400906S	9	60	6
4400908S	9	80	7
4400910S	9	100	7
4401002S	10	20	6
4401003S	10	30	6
4401004S	10	40	6
4401006S	10	60	7
4401008S	10	80	7
4401010S	10	100	7
4401202S	12	20	7
4401203S	12	30	7
4401204S	12	40	7
4401206S	12	60	7

135cm Shaft Length

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	CSI Fit F
4400304X	3	40	5
4400306X	3	60	5
4400308X	3	80	5
4400310X	3	100	5
4400312X	3	120	5
4400315X	3	150	5
4400322X	3	220	5
4400402X	4	20	5
4400404X	4	40	5
4400406X	4	60	5
4400408X	4	80	5
4400410X	4	100	5
4400412X	4	120	5
4400415X	4	150	5
4400422X	4	220	5
4400502X	5	20	5
4400504X	5	40	5
4400506X	5	60	5
4400508X	5	80	5
4400510X	5	100	5
4400512X	5	120	5
4400515X	5	150	5
4400522X	5	220	5
4400602X	6	20	5
4400604X	6	40	5
4400606X	6	60	5
4400608X	6	80	5
4400610X	6	100	5

135cm Shaft Length Continued

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	CSI Fit F
4400612X	6	120	5
4400615X	6	150	5
4400622X	6	220	5
4400702X	7	20	5
4400704X	7	40	5
4400706X	7	60	5
4400708X	7	80	5
4400710X	7	100	6
4400802X	8	20	5
4400804X	8	40	5
4400806X	8	60	5
4400808X	8	80	5
4400810X	8	100	6
4400902X	9	20	6
4400904X	9	40	6
4400906X	9	60	6
4400908X	9	80	7
4400910X	9	100	7
4401002X	10	20	6
4401004X	10	40	6
4401006X	10	60	7
4401008X	10	80	7
4401010X	10	100	7
4401202X	12	20	7
4401204X	12	40	7
4401206X	12	60	7

Advancing Endovascular Therapy

Cordis is making advancements in Endovascular therapy with the new Cordis POWERFLEX® Pro PTA Dilatation Catheter. Now featuring lengths up to 220 mm, we are offering a broad portfolio of specialized devices for treating complex and diffuse lesions.

PROFILE	POWER	PERFORMANCE
6 x 220 mm Through a 5F	RBP up to 18ATM	Short Shoulders for Accuracy
8 x 80 mm Through a 5F	Non-Compliant for Post Stent Dilatation	Rapid Inflation and Deflation



Important Information: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedures, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

Tel: +32 2 746 3000 Fax: +32 2 746 3405

© Cordis Europe, a division of Johnson & Johnson Medical N.V. - March 2012 2E-800-2000-2