

Achilles NC

0.018" system

Combination of strength and optimal control in routine and challenging cases

Broad range of sizes available

Confident to treat all lower extremity lesions

Non-Compliant OTW

Excellent for high pressure dilatation and lesion preparation

Seamless tip-to-shaft design

Provides smooth transmission and maximized pushability

Low tip profile

Facilitates lesion entry



Optimized non-compliant balloon material

High pressure balloon inflation up to 22atm

Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	10	15	20	30	40	60	80	100	150				
3.0	—	—	V	V	V	V	—	—	—	40	4	12	22
3.5	—	—	V	V	V	V	—	—	—	40	4	12	22
4.0	—	—	V	V	V	V	—	—	—	40	4	12	22
4.5	—	—	V	V	V	V	—	—	—	40	5	12	20
5.0	—	—	V	V	V	V	—	—	—	40	5	12	20
5.5	—	—	V	V	V	V	—	—	—	40	5	12	20
6.0	—	—	V	V	V	V	—	—	—	40	6	12	20
7.0	—	—	V	V	V	V	—	—	—	40	6	12	18
8.0	—	—	V	V	V	V	—	—	—	40	6	12	18
2.0	V	V	V	V	V	V	V	V	V	70	4	12	22
2.5	V	V	V	V	V	V	V	V	V	70	4	12	22
3.0	V	V	V	V	V	V	V	V	V	70	4	12	22
3.5	V	V	V	V	V	V	V	V	V	70	4	12	22
4.0	V	V	V	V	V	V	V	V	V	70	4	12	22
4.5	V	V	V	V	V	V	V	V	V	70	5	12	20
5.0	V	V	V	V	V	V	V	V	V	70	5	12	20
5.5	V	V	V	V	V	V	V	V	V	70	5	12	20
6.0	V	V	V	V	V	V	V	V	V	70	6	12	20
7.0	V	V	V	V	V	V	V	V	V	70	6	12	18
8.0	V	V	V	V	V	V	V	V	V	70	6	12	18
2.0	V	V	V	V	V	V	V	V	V	90	4	12	22
2.5	V	V	V	V	V	V	V	V	V	90	4	12	22
3.0	V	V	V	V	V	V	V	V	V	90	4	12	22
3.5	V	V	V	V	V	V	V	V	V	90	4	12	22
4.0	V	V	V	V	V	V	V	V	V	90	4	12	22
4.5	V	V	V	V	V	V	V	V	V	90	5	12	20
5.0	V	V	V	V	V	V	V	V	V	90	5	12	20
5.5	V	V	V	V	V	V	V	V	V	90	5	12	20
6.0	V	V	V	V	V	V	V	V	V	90	6	12	20
7.0	V	V	V	V	V	V	V	V	V	90	6	12	18
8.0	V	V	V	V	V	V	V	V	V	90	6	12	18
2.0	V	V	V	V	V	V	V	V	V	150	4	12	22
2.5	V	V	V	V	V	V	V	V	V	150	4	12	22
3.0	V	V	V	V	V	V	V	V	V	150	4	12	22
3.5	V	V	V	V	V	V	V	V	V	150	4	12	22
4.0	V	V	V	V	V	V	V	V	V	150	4	12	22
4.5	V	V	V	V	V	V	V	V	V	150	5	12	20
5.0	V	V	V	V	V	V	V	V	V	150	5	12	20
5.5	V	V	V	V	V	V	V	V	V	150	5	12	20
6.0	V	V	V	V	V	V	V	V	V	150	6	12	20
7.0	V	V	V	V	V	V	V	V	V	150	6	12	18
8.0	V	V	V	V	V	V	V	V	V	150	6	12	18

Notes: Sizes marked “V” means available, Sizes marked “—” have longer delivery time.



BrosMed Medical Co., Ltd.

Add: 15th Building, SMEs Venture Park SongShan Lake
 Hi-Tech Industrial Development Zone Dongguan
 523808, GuangDong, China
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
 E-mail: sales@brosmed.com www.brosmed.com

BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
 Tel: +31 (0) 850 140 901
 E-mail: cs@brosmed.com
 Web: www.brosmed.com

US FDA 510K cleared



2797
 GR4-6021/Rev.01