

Atropos sc

0.035" system

Combines strength and optimal control in routine and challenging cases

To treat fibrotic and calcified lesions

Broad size range:
Diameter: 3.0 - 10.0mm
Length: 20 - 200mm

Semi-Compliant OTW

Engineered for AV access

Excellent pushability


Seamless tip-to-shaft design for smooth push force transmission

Lowest tip profile 0.038"

Facilitates lesion crossing with ease

Hydrophilic coating

Latest hydrophilic coating technology to facilitate easy maneuvering through complex anatomy and lesions



Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	120	150	200				
3.0	V	V	V	V	—	—	—	—	—	40	5	6	14
4.0	V	V	V	V	—	—	—	—	—	40	5	6	14
5.0	V	V	V	V	—	—	—	—	—	40	5	6	14
6.0	V	V	V	V	—	—	—	—	—	40	6	6	14
7.0	V	V	V	V	—	—	—	—	—	40	6	6	14
8.0	V	V	V	V	—	—	—	—	—	40	6	6	14
3.0	V	V	V	V	V	V	V	V	V	70	5	6	14
4.0	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	70	6	6	14
7.0	V	V	V	V	V	V	V	V	V	70	6	6	14
8.0	V	V	V	V	V	V	V	V	—	70	6	6	14
9.0	V	V	V	V	V	V	V	—	—	70	7	6	14
10.0	V	V	V	V	V	—	—	—	—	70	7	6	14
3.0	V	V	V	V	V	V	V	V	V	90	5	6	14
4.0	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	90	6	6	14
7.0	V	V	V	V	V	V	V	V	V	90	6	6	14
8.0	V	V	V	V	V	V	V	V	—	90	6	6	14
9.0	V	V	V	V	V	V	V	—	—	90	7	6	14
10.0	V	V	V	V	V	—	—	—	—	90	7	6	14
3.0	V	V	V	V	V	V	V	V	V	150	5	6	14
4.0	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	150	6	6	14
7.0	V	V	V	V	V	V	V	V	V	150	6	6	14
8.0	V	V	V	V	V	V	V	V	—	150	6	6	14
9.0	V	V	V	V	V	V	V	—	—	150	7	6	14
10.0	V	V	V	V	V	—	—	—	—	150	7	6	14

Notes: Sizes marked “V” means available, Sizes marked “—” have longer delivery time.



BrosMed Medical Co., Ltd.

Add: 15th Building, SMEs Venture Park SongShan Lake
 Hi-Tech Industrial Development Zone Dongguan
 523808, GuangDong, China
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
 E-mail: sales@brosmed.com www.brosmed.com

BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
 Tel: +31 (0) 850 140 901
 E-mail: cs@brosmed.com
 Web: www.brosmed.com

US FDA 510K cleared



2797
 GRA-6022/Rev.01