



# Hemostasis Valves & Accessories



Minimizing blood loss and optimizing control of interventional devices.



## AccessPLUS™

The AccessPLUS large bore hemostasis valve is a second generation device that maximizes the “sweet spot” or optimal sealing surface for blood control.

### SPECIFICATIONS

Inner lumen: .120” (3.1mm)  
Body material: Polycarbonate



## DoublePlay™

Two large bore hemostasis ports allow easier manipulation and control of interventional devices during “kissing balloon” techniques.

### SPECIFICATIONS

Inner lumen: .120” (3.1mm)  
Body material: Polycarbonate



## Passage®

The Passage represents the first of Merit’s standard line of hemostasis valves. It accommodates devices up to 7.5 French.

### SPECIFICATIONS

Inner lumen: .097” (2.5mm)  
Body material: Polycarbonate



## Honor®

This is the latest in Merit’s hemostasis valve technology. The easy push-pull mechanism of the Honor Hemostasis Valve allows single-handed manipulation to open and close the valve while introducing devices into the vasculature. The hemostasis seal minimizes blood loss while giving the clinician the ability to manipulate interventional devices.

### SPECIFICATIONS

Inner lumen: .094” (2.4mm)  
Body material: Polycarbonate



## MBA™ - Merit’s Bleedback Alternative™

A premier hemostasis valve which features two seals. A conventional “touhy” silicone hemostasis valve provides security that the interventional device will stay in place. The secondary slit seal delivers a virtually bloodless procedure.

### SPECIFICATIONS

Inner lumen: .118” (3mm)  
Body material: Polycarbonate



## FLO30™

Merit Medical’s FLO30™ hemostasis valve is a single-use, proprietary valved adapter which provides a leak-proof seal during interventional and diagnostic procedures.

### SPECIFICATIONS

Inner lumen: .118” (3mm)  
Body material: Polycarbonate

## Merit Accessory Pak

In response to the growing demand for minimized blood loss, reduced exposure to blood, and easy to use devices, Merit is the first to make available a full line of hemostasis valves with unequalled options and sizes necessary to meet your procedural demands and professional preferences.



**Hemostasis valves and accessories are also available with any Merit inflation device. Kits are available in standard catalog numbers as well as custom kit configurations**

## BasixTOUCH™

The BasixTOUCH™ Inflation device is an innovative device with an ergonomic handle that enables one-handed preparation and priming, 35 atm pressure capability and 30 mL fluid capacity. Features a quick-release handle for easy release, even from high-pressure inflations, accelerating deflation time. Reaches high pressures rapidly—with less than three rotations of the handle. Comes with high pressure stopcock.



## Blue Diamond®

The 30 atm Blue Diamond inflation device combines the accuracy of a digital device with the simplicity of analogue systems. Features include one-button functionality for simplicity, a large dial face with oversized digits to enhance visibility, and last inflation marks to recall the highest pressure of the previous inflation. Comes with 13" (33 cm) extension tubing and 3-way stopcock. Packaged 5 units per box, 4 boxes per case.



## BasixCOMPAK™

The BasixCOMPAK is a 30 atm analogue inflation device that features an angled gauge for better viewing of pressure readings and a luminescent dial for low-light procedures.



## BasixTOUCH40™

The BasixTOUCH40 has a 40 ATM pressure capacity allows inflation of high pressure interventional balloons. 30 mL syringe barrel accommodates additional fluid for large balloons.



## ORDERING INFORMATION



### Large Bore 9 French Inner Lumen

Found in the MBA, Access-9, AccessPLUS, DoublePlay, and FLO30 designs.

### Airless Rotator

A proprietary airless rotator minimizes the possibility of trapped air.

Catalog Number	Hemostasis valve	Accessories		
		Metal wire introducer	Plastic wire introducer	Torque device

#### Standalone Hemostasis Valves

MAP150	Access plus, large bore			
MAP220	Double Play, large bore			
MAP300	Honor, small bore			
MAP350	Honor, small bore with bonded 20cm tubing and a stopcock			
MAP400	MBA, large bore			
MAP450	MBA, large bore with bonded 25.4cm high pressure tubing			
K12-00040	Passage, small bore with bonded 15.2cm high pressure tubing			
FLO30	FLO30, large bore			

#### Hemostasis Valve MAP Kits

MAP101	Passage, small bore	X		
MAP102	Passage, small bore	X		X
MAP111	Access 9, large bore	X		
MAP112	Access 9, large bore	X		X
MAP113	Access 9, large bore		X	
MAP151	Access plus, large bore	X		
MAP152	Access plus, large bore	X		X
MAP153	Access plus, large bore		X	
MAP154	Access plus, large bore		X	X
MAP221	Double Play, large bore	X		
MAP222	Double Play, large bore	X		X
MAP301	Honor, small bore	X		
MAP302	Honor, small bore	X		X
MAP351	Honor, small bore with bonded 20cm tubing and a stopcock	X		
MAP352	Honor, small bore with bonded 20cm tubing and a stopcock	X		X
MAP353	Honor, small bore with bonded 50cm tubing and a stopcock	X		X
MAP403	MBA, large bore	X		X
MAP453	MBA, with bonded 25.4cm high pressure tubing	X		X

#### Accessories

MAP500	Torque Device 0.018", salmon			X
MAP501	Torque Device 0.018", green			X
MAP550	Metal Wire Insertion Tool	X		
MAP600	Plastic Wire Insertion Tool		X	

Packaged 25 units per box.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

merit.com

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Korea  
+8210 2205 0707

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 91 1238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115