

# Polux SC

## 0.014" system

Excellent performance in routine and challenging cases

Designed to treat resistant lesions

For a wide range of indications

## Semi-Compliant OTW

### Designed to treat BTK lesions

### Optimized crossing profile

Best in class crossing capability for diffused and challenging lesions

### Hydrophilic coating


Ensures smooth entry to tight lesions

### Kink resistant shaft

Smooth transmission through tortuous vasculature

### Low entry profile 0.017"

For reaching distal vasculature



# Specifications

Balloon Diameter (mm)	Balloon Length (mm)												Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	5	10	15	20	30	40	60	80	100	120	150	200				
1.5	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	70	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	70	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	90	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	90	5	6	14
1.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
2.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
3.5	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.0	V	V	V	V	V	V	V	V	V	V	V	V	150	4	6	14
4.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
5.5	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	V	V	V	150	5	6	14

Notes: Sizes marked “V” means available, Sizes marked “—” have longer delivery time.



## BrosMed Medical Co., Ltd.

Add: 15th Building, SMEs Venture Park SongShan Lake  
 Hi-Tech Industrial Development Zone Dongguan  
 523808, Guangdong, China  
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016  
 E-mail: sales@brosmed.com www.brosmed.com

## BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands  
 Tel: +31 (0) 850 140 901  
 E-mail: cs@brosmed.com  
 Web: www.brosmed.com

US FDA 510K cleared



GRA-6023/Rev.01